




Product Special

Rib Spreader Cooley Mitral Rippenspreizer Cooley Mitral



Product information Produktinformationen

Rippenspreizer Cooley Mitral

- komplett mit Zubehör nach Cosgrove
- verfügbar in 2 Größen (XL und XXL)
- Edelstahl

Rib Spreader Cooley Mitral

- complete with Cosgrove mounting bars
- available in 2 sizes (XL and XXL)
- stainless steel



Unsere neue Produkt-Welt!
Our new Produc-World!

made in GERMANY

Änderungen der technischen Angaben und Maße liegen in unserem Ermessen. Modelländerungen vorbehalten. Ohne unsere ausdrückliche schriftliche Genehmigung darf dieses Printmedium weder im Gesamten noch auszugsweise kopiert oder nachgedruckt werden.

Changes in technical specifications and dimensions are at our discretion. Subject to modification. Without our explicit written approval this print media may not be copied or reprinted as a whole or in extracts.

COOLEY MITRAL SET XXL

1-10-708-40	according to Cosgrove, consisting of:			
	Type	spreading width	arm length	Qty
1-10-708-42	Cooley Mitral Frame, XXL	211 mm	218 mm	1
1-10-708-45	Side arm bar, blade retractor attachment, diam: 6 mm			1
1-10-708-50	Main clamp blade retractor attachment, diam: 6 mm			1
1-10-729-30	Universal clamp, clamps both sides			3
	Type	width	depth	Qty
1-10-709-05	Cooley blade small, curved	28 mm	35 mm	1
1-10-709-15	Cooley blade medium, curved	38 mm	35 mm	1
1-10-709-20	Cooley blade flat	42 mm	22 mm	3

OPTIONS


	Type	spreading width	arm length	Qty
1-10-708-41	Cooley Mitral Frame XL	187 mm	166 mm	1
	Type	width	depth	Qty
1-10-708-60*	Cooley atrial retractor blade, small	25 mm	35 mm	1
1-10-708-65*	Cooley atrial retractor blade	40 mm	35 mm	1
1-10-708-70*	Cooley atrial retractor blade, large	50 mm	35 mm	1
1-10-708-75*	Cooley atrial retractor blade, tricuspid	50 mm	20 mm	1
1-10-729-45	Atrial Retractor, large	50 mm	35 mm	1
1-10-729-48	Atrial Retractor, large, low profile	50 mm	20 mm	1
1-10-729-50	Ribbon blade retractor	25 mm	76 mm	1
1-10-729-52	Ribbon blade retractor	19 mm	76 mm	1
1-10-729-55	Ribbon blade retractor	12,7 mm	76 mm	1
1-10-729-35	Cooley small blade	25 mm	35 mm	1
1-10-709-10	Cooley blade medium, slightly curved	38 mm	48 mm	1
1-10-729-40	Cooley medium blade	38 mm	35 mm	1
1-10-709-25*	Cooley blade small, curved	28 mm	35 mm	1
1-10-709-30*	Cooley blade medium, slightly curved	38 mm	48 mm	1
1-10-709-35*	Cooley blade medium, curved	38 mm	35 mm	1
1-10-709-40*	Cooley blade flat	42 mm	22 mm	1

*malleable



Product Special

Rib Spreader Sellors Rippenspreizer Sellors



Product information Produktinformationen

Rippenspreizer Sellors

- gebogene Arme
- schlüssellose Adaption der Valven an den Rahmen
- Edelstahl

Rib Spreader Sellors

- curved arms
- keyless adaption of blades to frame
- stainless steel



Unsere neue Produkt-Welt!
Our new Product-World!

made in **GERMANY**

Änderungen der technischen Angaben und Maße liegen in unserem Ermessen. Modelländerungen vorbehalten. Ohne unsere ausdrückliche schriftliche Genehmigung darf dieses Printmedium weder im Gesamten noch auszugsweise kopiert oder nachgedruckt werden.

Changes in technical specifications and dimensions are at our discretion. Subject to modification. Without our explicit written approval this print media may not be copied or reprinted as a whole or in extracts.

SELLORS RIB SPREADER COMPLETE

1-10-739-05	consisting of:			
	Type	spreading width	arm length	Qty
1-10-739-10	Sellors Frame	min. 197 mm max. 226 mm	140 mm	1

	Type	width	depth	Qty
1-10-739-25	1 pair of blades	50 mm	40 mm	1
1-10-739-30	1 pair of blades	50 mm	50 mm	1
1-10-739-35	1 pair of blades	50 mm	70 mm	1

OPTIONS

	Type	width	depth	Qty
1-10-739-50	1 pair of blades	75 mm	30 mm	1
1-10-739-55	1 pair of blades	75 mm	45 mm	1
1-10-739-60	1 pair of blades	120 mm	30 mm	1
1-10-739-65	1 pair of blades	75 mm	78 mm	1